

Plastic Tube Container BlockPack

Square telescopic packaging tubes with twist lock closure feature and ratchet length adjustment

Material

→ Polyethylene (PE)

Product Range

→ PE natural

→ outer parts and inner parts are packed separately

Production environment

→ numerous sizes available in cleanroom production



Product features BlockPack

- highly secure bayonet-type twist lock closure, requires minimal force in opening and closing the pack
- contoured tube ends and lateral support grooves are integrated together to provide a cushioning effect
- with ratchet length adjustment
- ergonomic product design offer shape stability
- speedy handling
- safe transportation
- excellent for stacking and storage
- durable and reusable
- can be used as primary or secondary packaging

The proof of the suitability of the packaging in combination with the medical device to be packed must be provided by the manufacturer.










Product Range: BlockPack








Production Environment

 Cleanroom production  Industrial standard / Cleanroom production upon request

Article	Inner diam.	Usable length adjustable	Packaging Unit	Production Environment
BK 8 050	8mm	50 - 80mm	3000	
BK 8 080	8mm	80 - 120mm	3000	
BK 8 120	8mm	120 - 200mm	1500	
BK 8 200	8mm	200 - 350mm	2000	
BK 10 050	10mm	50 - 80mm	2100	
BK 10 080	10mm	80 - 120mm	2000	
BK 10 120	10mm	120 - 200mm	1170	
BK 10 200	10mm	200 - 350mm	1100	
BK 12 050	12mm	50 - 80mm	2200	
BK 12 080	12mm	80 - 120mm	1550	
BK 12 120	12mm	120 - 200mm	800	
BK 12 200	12mm	200 - 350mm	800	
BK 12 250	12mm	250 - 460mm	700	
BK 14 050	14mm	50 - 80mm	1700	
BK 14 080	14mm	80 - 120mm	1100	
BK 14 120	14mm	120 - 200mm	800	
BK 14 200	14mm	200 - 350mm	600	
BK 16 050	16mm	50 - 80mm	1600	
BK 16 080	16mm	80 - 120mm	850	
BK 16 120	16mm	120 - 200mm	1200	
BK 16 200	16mm	200 - 350mm	1200	
BK 18 050	18mm	50 - 80mm	1300	
BK 18 080	18mm	80 - 120mm	850	
BK 18 120	18mm	120 - 200mm	1050	
BK 18 200	18mm	200 - 350mm	400	
BK 20 050	20mm	50 - 80mm	1200	
BK 20 080	20mm	80 - 120mm	1300	
BK 20 120	20mm	120 - 200mm	900	
BK 20 200	20mm	200 - 350mm	800	
BK 20 350	20mm	350 - 620mm	500	
BK 22 050	22mm	50 - 80mm	1000	
BK 22 080	22mm	80 - 120mm	1200	
BK 22 120	22mm	120 - 200mm	800	
BK 22 200	22mm	200 - 350mm	600	
BK 24 050	24mm	50 - 80mm	750	
BK 24 080	24mm	80 - 120mm	1000	
BK 24 120	24mm	120 - 200mm	700	

Article	Inner diam.	Usable length adjustable	Packaging Unit	Production Environment
BK 24 200	24mm	200 - 350mm	550	
BK 26 050	26mm	50 - 60mm	700	
BK 26 060	26mm	60 - 80mm	700	
BK 26 080	26mm	80 - 120mm	850	
BK 26 120	26mm	120 - 200mm	600	
BK 26 200	26mm	200 - 350mm	450	
BK 26 350	26mm	350 - 620mm	288	
BK 32 050	32mm	50 - 70mm	1200	
BK 32 060	32mm	60 - 80mm	800	
BK 32 080	32mm	80 - 120mm	600	
BK 32 120	32mm	120 - 200mm	400	
BK 32 200	32mm	200 - 350mm	360	
BK 32 350	32mm	350 - 620mm	140	
BK 35 050	35mm	50 - 70mm	900	
BK 35 060	35mm	60 - 80mm	700	
BK 35 080	35mm	80 - 120mm	500	
BK 35 120	35mm	120 - 200mm	300	
BK 35 200	35mm	200 - 350mm	300	
BK 35 350	35mm	350 - 620mm	159	
BK 40 080	40mm	80 - 120mm	450	
BK 40 120	40mm	120 - 200mm	300	
BK 40 200	40mm	200 - 350mm	240	
BK 40 350	40mm	350 - 620mm	120	
BK 45 060	45mm	60 - 80mm	450	
BK 45 080	45mm	80 - 120mm	350	
BK 45 120	45mm	120 - 200mm	275	
BK 45 200	45mm	200 - 350mm	160	
BK 45 200	45mm	200 - 350mm	1000	
BK 45 350	45mm	350 - 620mm	90	
BK 55 060	55mm	60 - 80mm	250	
BK 55 080	55mm	80 - 120mm	200	
BK 55 120	55mm	120 - 200mm	200	
BK 55 200	55mm	200 - 350mm	100	
BK 55 350	55mm	350 - 620mm	70	
BK 65 080	65mm	80 - 120mm	210	
BK 65 120	65mm	120 - 200mm	120	
BK 65 200	65mm	200 - 350mm	90	

Article	Inner diam.	Usable length adjustable	Packaging Unit	Production Environment
BK 65 350	65mm	350 - 620mm	48	
BK 75 080	75mm	80 - 120mm	168	
BK 75 120	75mm	120 - 200mm	96	
BK 75 200	75mm	200 - 350mm	58	
BK 85 080	85mm	80 - 120mm	90	
BK 85 120	85mm	120 - 200mm	72	
BK 85 200	85mm	200 - 350mm	48	
BK 85 350	85mm	350 - 620mm	30	
BK 85 550	85mm	550 - 1000mm	20	

Article	Inner diam.	Usable length adjustable	Packaging Unit	Production Environment
BK 95 080	95mm	80 - 120mm	80	
BK 95 120	95mm	120 - 200mm	60	
BK 95 200	95mm	200 - 350mm	30	
BK 95 350	95mm	350 - 620mm	20	
BK 105 080	105mm	80 - 120mm	80	
BK 105 120	105mm	120 - 200mm	48	
BK 105 200	105mm	200 - 350mm	30	
BK 105 350	105mm	350 - 620mm	20	